

## MEGNETIC STIRER WITH HOT PLATE



s.no	Particulars	compliance
1	Material of construction	Body made of mild steel duly powder coated. Top plate of stainless steel.
2	Temp Range	Ambient to 300c.
3	Stirrer RPM	1400 controlled by steeples regulator
4	Temp. Controller	Digital microprocessor based with thermocouple which can be inserted in the material.
5	RPM indicator	Digital (Optional)
6	Operating Voltage	220-240 Volts Ac
7	Load	1kw max.
8	Capacity	5 lts. Of water
9	Stirrer bar	Provide of size 25mm & 38mm