


Solvent Filtration Kit

Solvent Filtrations Assembly Glass / SS (upper part) with 1 lit. Bottom Flask Glass Size 47mm

It used to filter the Solvent (mobile phase) for HPLC & other analytical application. Upper part is of Glass or SS as per requirement.



Model: SFK-GC



Model: SFK-SS




Vacuum Pump Model PCI-15
Max. Flow - 15 LPM,
Vacuum - 22" Hg (554 mm)
Max. Pressure - 25 Psig, 1/20 HP

Solvent Filtration in Mobile Phase Bottle



Solvent & Sample Filtration Kit



Nylon 66 Membrane Filter
Size: \varnothing 13 x 0.2 μ (Porosity)
Other sizes & porosity also available

\varnothing 13mm S. S. Filter Holder /
 \varnothing 13mm P. P. Filter Holder


5ml Glass Hypodermic Syringe (Indian) /
5ml Gas Tight Syringe (Indian)

Rheodyne Needle

6 Branch
Model: STV-6B

3 Branch
Model: STV-3B

Multi-Fold Solvent Filtration System



Model: PCI 75-F

Model: SFS-MF

Salient Features:

- Manifold with 3 ball valves (ON/OFF)
- Single Mobile Phase in triple quantity
- Three different mobile phases simultaneously (choice 1, 2 or 3)
- Single Vacuum Pump PCI-75S to operate 3 Solvent Filtration
- Ease of operation, maintenance free
- Improves efficiency in heavy load situations
- Innovative design by PCI Analytics only
- No Liquid contact with Manifold & hence no contamination
- Vacuum Flask 1 Ltr., 2 Ltr. capacity