



An ISO 9001:2008 Certified Company

CREATING STANDARDS  
WITH COMMITMENT TO QUALITY

## LABORATORY OVEN (UNIVERSAL MEMMERT TYPE)

Double walled construction. Inner chamber made of anodized Aluminum / S.Steel. Outer body made of mild steel is powder coated/painted. Temperature thermostatically controlled by hydraulic thermostat from 50°C-250°C. Elements are at the bottom and sides. Built in 'L' shape thermometer. The gap between the inner & outer chamber is filled with high grade glass wool to avoid the thermal losses. Perforated adjustable trays. Two pilot lamp on/off switch. To work on 220/230 V.AC.

SHI-102



Model No.	Inner Chamber (W x D x H)	Volume (Liters)	Load (KW)	Shelves (No's)
SHI-101a	300 x 300 x 300 mm	28	1.5	2
SHI-101B	355 x 355 x 355 mm	43	1.5	2
SHI-101C	455 x 455 x 455 mm	90	2.0	2
SHI-101D	455 x 455 x 605 mm	105	2.2	3
SHI-101E	605 x 605 x 605 mm	215	2.5	3
SHI-101F	605 x 455 x 910 mm	240	2.5	3
SHI-101G	605 x 605 x 910 mm	325	3.0	3

Other Custom Sizes are also Available.....

--	--	--	--	--

📞 9891673394, 9716963092

✉ banditech07@gmail.com | info@banditechnology.com 🌐 www.indiamart.com/banditechnology