



## AUTOCALVE VERTICAL CE



Outer body made of thick mild steel finished oven baked enamel. Inner chamber will be thick stainless steel. The chamber is absolutely leak-proof and it can be operated at any selective point in between 5 to 20 pound/Sq. inch. The instrument is equipped with pressure gauge. Pressure release valve, safety valve, water level indicator and perforated S.S. basket complete with cord and plug. The lid is made of thick mild steel sheet duly lined with stainless steel with pedal lifting arrangements and radial locking devices to work on 220/230V AC single phase supply

Inner Chamber Model	Load capacity	Outer M.S.	Outer body SS	GMP
---------------------	---------------	------------	---------------	-----

Dia x Depth (Extra Cost)

a) 250 x 450 mm	1.5 kw	22
b) 300 x 500 mm	2.0 kw	40
c) 350 x 550 mm	3.0 kw	50
d) 450 x 600 mm	4.0 kw	98
e) 550 x 750 mm	5.0 kw	152

Optional Accessories

1. Automatic low water level cut off device
2. Pressure control switch
3. Automatic 0-60 minute timer
4. Digital Temperature indicator
5. Automatic control system for temp. Pressure

**Other Custom Sizes are also Available.....**

☎ 9891673394, 9716963092

✉ banditech07@gmail.com | info@banditechnology.com 🌐 www.indiamart.com/banditechnology