


SHI-174



SHI-173



### **SHI-173 : SIEVE SHAKER (GRYTORY MOTORIZED)**

This can accommodate up to 7 sieves of 6" to 8" dia and is driver by ¼ H.P. Electric Motor through reduction gear. Every part of the sieves is utilized by the constant upward and downward movement. The whole gear mechanism runs in oil bath supplied without test sieves suitable for operation on 220 Volt AC.

### **OPTIONAL**

- do- fitted with imported timer of 0-60 minutes
- do- fitted with crompton / equillent motor 1/4 H.P.

### **SHI-174 : SIEVE SHAKER (ROTOP MOTORIZED)**

This is heavy machine driven by ¼ H.P. Motor. The circular end to end motion of 250 minimum to give by mechanism which is fully immersed in oil. The tapping motion of 150 is given by hemmer from the top. The shaker takes 6" to 8" sieves. It is recommended of ISS or BSS test sieves, the top, the shaker takes 6" to 8 sieves. It is recommended of ISS or BSS test sieves.

### **OPTIONAL:**

- do- fitted with imported timer 0-60 minutes
- do- fitted with crompton /equillent motor ¼ H.P. at extra price

## TEST SIEVES



### TEST SIEVES FINE SERIES IN 200MM DIA BRASS FRAME

BSS	ASTM	ISI	MICRONS	MM	
4	5	--		4.75	
5	6	340	1812	3.35	
6	7	280	2411	2.80	
7	8	240	2207	2.36	
8	10	200	2057	2.00	
10	12	170	1698	1.70	
12	14	140	1005	1.40	
14	16	120	1200	1.18	
16	18	100	1000	1.00	
18	20	85	850	0.850	
22	25	70	710	0.710	
25	30	60	600	0.600	
30	35	50	500	0.500	

# BANDI TECHNOLOGY

Note : Contact for Labrotary Equipment

Mob: +91-9891673394  
+91-9716963092  
e-mail : banditech07@gmail.com

36	40	40	425	0.425	
44	45	35	355	0.355	
52	50	30	300	0.300	
60	60	25	250	0.250	
72	70	20	212	0.212	
85	80	18	180	0.180	
100	100	15	150	0.150	
120	120	12	125	0.125	
150	140	10	106	0.106	
170	170	9	90	0.90	
200	200	8	75	0.75	
240	230	7	63	0.63	
300	270	6	53	0.53	
325	300	5	45	0.45	
350	325	4	44	0.044	
400	380	3	40	0.040	
500	---	--	25	0.025	

## **TEST TIEVES OF COURSE SERIES:-**

Made of perforated plate sheet with square holes. Available in 200 /300/450mm, dia.

Aperture size:-

125,160,100,90,80,75,63,53,50,45,40,37.5,31.5,25,22.4,19,16,13.2,12.2,11.2,10.2,10,9.5,8,6,7,6.3,  
5.6,4.75,2.8mm

S-494, 3<sup>rd</sup> Floor, School Block, Shakarpur, Delhi-110092

[www.banditechnology.com](http://www.banditechnology.com)