





## BALL MILL

 An ISO 9001:2015 Certified Company  
ISO 13485:2016

### SHI-176:LAB. WILLYMILL (ARTHURH. THOMASTYPE)

Willy mill is fitted on a substantial cast iron base. Fitted with suitable motor supplied with S.S. HOPPER with cover, three delivery sieves tops of 20, 40 and 80 mesh with plastic receiver, brush,. A removable Plexi plate permits observation during operation. This Unit comes with two stationary steel blades and rotary with steel edge provided in the cast iron chamber Suitable to operate on 220 Volts, Single phase, 50 Hz cycles A.C. Supply



**Applications**

- Grinding of small samples
  - Poultry Farms
  - Tobacco Industry
- Biological Research Lab
  - Agriculture
  - Botanical Research Lab

### Technical Specifications:-

S.NO	PARTICUALRS	COMPLIANCE	
1	MOC	Body is made of mild steel powder coated &fitted on heavy Base .	
2	Contact Parts	Made of stainless steel 304 Grade.	
3	Speed Regulator	700 RPM	
4	Capacity chamber size	Grinding Chamber size	Load
		40x25mm	1/4 hp
		65x25mm	1/4 hp
		100x25mm	1/2 hp
200x7500	1hp		
5	Delivery Output size	20,40,60 mesh size	
7	Voltage	220-240 volts 50 HZ	
8	Spares	Brush, sieves of 20,40,60 mesh , spare jars.	

📞 9891673394, 9716963092

✉ banditech07@gmail.com | info@banditechnology.com 🌐 www.indiamart.com/banditechnology